



# OPEN CHANNEL DRAINAGE (BLAKE TYPE) CODE 125




This product has been designed for clearing bloody, serous or purulent fluid or air collections from body compartments as a result of post-surgical waste (abdominal, pelvic, retroperitoneal, cardiovascular, chest and plastic surgery).

- Circular section of drainage tube (80 cm long).
- Long-term indwelling.
- Larger area of contact with secretions.
- Rapid and efficient absorption.
- Improved patient comfort.



Product  
**100%**  
Silicone



100% clear silicone with full-length radiopaque line.

Longitudinal absorption along the complete inserted tube.



Step connector.

Code	Description	Channelled section length (cm)	Total length (cm)
125-10	10.0 FR open channel drainage	30,00	210,00
125-15	15.0 FR open channel drainage	30,00	210,00
125-19	19.0 FR open channel drainage	30,00	210,00