



# FLAT DRAINAGE (JACKSON PRATT TYPE) CODE127



**SILMAG**  
Medicine evolves,  
so do we.

This product has been designed for clearing bloody, serous or purulent fluid or air collections from body compartments as a result of post-surgical waste (abdominal, pelvic, retroperitoneal, cardiovascular, chest and plastic surgery).

- Greater absorption capacity.
- Flat tube with radiopaque line.
- Long-term indwelling.
- Improved patient comfort.
- Easy to handle.
- Lower risk of flow obstruction due to kinking and/or crushing.



Product  
**100%**  
Silicone



*Multi-perforated flat section.*



*Grooved interior. Improves fluid drainage at any point along the flat segment.*

Code	Width (mm)	Long fenestrated (cm)	Long total (cm)
127-07	7,00	16,00	77,00
127-10	10,00	21,20	112,00
127-13	13,00	24,70	114,50