




EXTERNAL VENTRICULAR DRAINAGE

CODE 130



Appropriate for cerebrospinal fluid (CSF) drainage and for cases where of intracranial pressure (ICP) monitoring is required; for ventricle decompression in head trauma patients, in patients with cranioencephalic trauma and in stroke patients, and to reduce intracranial pressure (ICP).

- Luer lock connector.
- Guidewire - stainless steel wire.
- The closed system minimizes the risk of infection.
- Easy to regulate drainage..



Code	Size (FR)	Outer Ø (mm)	Inner Ø (mm)	Length (cm)
130-09	9,0	3,00	1,75	30,00