



# CHEST DRAINAGE TUBE

CODE150



This product has been designed for clearing bloody, serous or purulent fluid or air collections from body compartments as a result of post-surgical waste (chest surgery or trauma).

- Biocompatible.
- Side holes at the distal end.
- Improved patient comfort.
- Lower risk of unnecessary procedures.



*Radiopaque line.*



*Open, rounded end.*

Code	Size (FR)	Outer Ø (mm)	Inner Ø (mm)	Length (cm)	Holes
150-28	28,0	9,30	6,40	150,00	1
150-40	40,0	13,30	9,50	180,00	1
150-40-01	40,0	13,30	9,50	180,00	4