



# MUCUS SUCTION TUBE

## CODE 210



**SILMAG**  
Medicine evolves,  
so do we.

This product has been designed for aspirating secretions with an open system in patients with reduced cough reflex who present with secretions in the upper respiratory tract, intubated patients and patients with tracheostomy.

- Does not cause friction injury.
- Easy to handle.



100% silicone,  
45 cm long tube.



Open, rounded end,  
with 4 side perforations.



Regulator  
of suction.

Version with cap  
only for code 210-06,  
The rest of the versions  
come without a cap

Code	Size (FR)	Outer Ø (mm)	Inner Ø (mm)
210-06	6,0	2,00	1,00
210-08	8,0	2,70	1,50
210-10	10,0	3,30	2,00
210-12	12,0	4,00	2,50
210-14	14,0	4,70	3,00
210-16	16,0	5,30	3,50
210-18	18,0	6,00	4,00