



CPAP SYSTEM

CODE 233




This product has been designed to administer a mixture of air/oxygen with positive pressure in spontaneously breathing neonatal patients.

Kit:

- (A) 2 corrugated tubes of 140 cm.
- (B) Pressure monitoring line (PEEP).
- (C) 2 adapters for humidification from 22 mm to 10 mm.
- (D) 1 luer lock lid for closing.
- (E) 1 double slip luer adapter.
- (F) 1 standard 15 mm connector.
- Fixation system for the cannula:
- (G) 1 thread cap (cotton).
- (H) 2 insurance with double locking.
- (I) 2 fastening straps
- (J) Hydrocolloid and velcro.

- Does not cause pressure injury to the nasal septum.
- Adaptable to bubble CPAP or ventilator.
- Replacement parts help reduce costs.
- Replacement kit: cannula, connectors and cap.
- Corrugated tubes (allow for mobility without risk of kinking).



100% clear silicone cannula.
Adaptable to nostrils.

Code	Size (FR)	Cannula Ø (mm)	Distance between cannulas (mm)	Cannula length (mm)	Recommendation based on patient weight	Color bonnet
233-00	8,0	2,70	2,70	9,00	Menos de 1000 g	Yellow
233-01	10,0	3,30	2,70	9,00	De 1000 a 1500 g	Green
233-02	12,0	4,00	4,00	10,00	De 1500 a 2000 g	Green
233-03	15,0	5,00	4,50	12,50	De 2000 a 3000 g	Brown