




# DOUBLE PIGTAIL URETERAL STENT CODE 405



This product has been designed to drain urine. In transient obstructive processes, in patients undergoing surgery ureteral lithotripsy, pathologies obstructive oncological, traumatic stenosis or surgical, congenital malformations ureteral, kidney transplant.

- Greater comfort and decreased pain due to material quality.
- Does not cause friction injury.
- Multiperforated.
- Radiopacity helps visualization and proper location.
- Lower risk of obstruction due to accumulated solutes.



- Kit:
- (A) 1 Polyurethane catheter.
  - (B) 1 clamp.
  - (C) 1 Teflon-coated wire.
  - (D) 1 Straight pusher.

Code	Size (FR)	Length (cm)
405-47-20	4,7	20,00
405-47-24	4,7	24,00
405-06-26	6,0	26,00
405-06-28	6,0	28,00
405-07-26	7,0	26,00
405-07-28	7,0	28,00

Open/closed end

Code	Size (FR)	Length (cm)
405-47-20-A	4,7	20,00
405-47-24-A	4,7	24,00
405-06-26-A	6,0	26,00
405-06-28-A	6,0	28,00
405-07-26-A	7,0	26,00
405-07-28-A	7,0	28,00

Open/open end

All named above, codes are available in hard and soft versions.