




CAVAL EPICUTANEOUS CATHETER CODE 540-541



This product has been designed for the administration of parenteral nutrition, hyperosmolar solutions, vasoactive drug infusion, intravenous solutions, in neonatal, infant and paediatric patients requiring prolonged parenteral treatment, patients with necrotizing enterocolitis (NEC), patients hospitalized in intensive care units, surgery patients requiring prolonged fasting, patients with metabolic alterations.

- Centimetric marking.
- At the end of the extension has a connection fitting that allows the catheter to be held simply, facilitating removal of the butterfly introducer.
- Ensures long-term indwelling and maximal X-ray visualization.



Kit:

- (A) 1 100% radiopaque silicone catheter.
- (B) 1 disposable tape,
- (C) 1 butterfly type introducer.

Caval epicutaneous catheter - FR 2.0			
Code	Length (cm)	Inner Ø x Outer Ø (mm)	* Flow rate (ml/min)
540-15	15,00	0,30 X 0,60	5,00
540-30	30,00	0,30 X 0,60	3,00

*Flow rate under 1 bar pressure.

Caval Epicutaneous Catheter - Premature - FR 1,8			
Code	Length (cm)	Ø Interior x Ø Exterior (mm)	* Flow rate (ml/min)
541-15	15,00	0,28 X 0,57	3,50
541-30	30,00	0,28 X 0,57	2,00

*Flow rate under 1 bar pressure.