



## DOUBLE LUMEN UMBILICAL CATHETER CODE 544-D



**SILMAG**  
Medicine evolves,  
so do we.

This catheter has been designed for placement in the umbilical vein of critically ill neonates hospitalized in neonatal intensive care units. Used for incompatible solutions or blood sample collection and in patients with blood type incompatibility requiring an exchange transfusion.

- Centimetre marks from 4 cm to 24 cm.
- Rounded end reduces the risk of vessel injury.
- Easy to handle and slide.
- Marking ensures positioning and maintenance during treatment.



*Full-length double lumen catheter for infusion of incompatible solutions.*



*Features a proximal end Luer lock.*



*100% radiopaque polyurethane and echogenic catheter.*

Code	Size (FR)	Outer Ø (mm)	Length (cm)
544-35-D	3,5	1,30	38
544-50-D	5,0	1,75	38