



PICC PERCUTANEOUS CATHETER WITH PEEL OFF INTRODUCER (DOUBLE LUMEN) CODE 650-D




This product has been designed for the administration of parenteral nutrition, hyperosmolar solutions, vasoactive drug infusion, intravenous solutions, in neonatal, infant and paediatric patients requiring prolonged parenteral treatment, patients with necrotizing enterocolitis (NEC), patients hospitalized in intensive care units, surgery patients requiring prolonged fasting, patients with metabolic alterations.

- Longer durability and lower catheter replacement rate.
- Lower cost per hospital day.
- Inner diameter ensures the right volume of fluid required for the patient.
- Enables administration of hyperosmolar solutions.
- Peel off introducer and side eye needle.
- Clear sheath for verification of venous blood return.



100% radiopaque silicone and echogenic catheter. 1 cm marks.



Double lumen



Double lumen facilitates administration of incompatible solutions.

Kit:

- (A) 2 needle-free connectors,
- (B) 1 disposable measuring tape,
- (C) 1 sterile silicone haemostatic tourniquet,
- (D) 1 metal forceps,
- (E) 1 peel off introducer,
- (F) 2 extensions.



Code	Size (FR)	Outer Ø (mm)	Length (cm)	Flow rate (ml/min)*
650-03-50-D	3,0	1,00	50	2x2,0
650-04-65-D	4,0	1,40	65	2x6,0
650-05-65-D	5,0	1,70	65	-

*Measurement made at 1 bar of pressure.