



PerfuPort
HUBER NEEDLE
CODE 1501



PerfuPort is a huber needle designed for fluid administration (perfusion) In patients with implantable catheters with titanium port or polysulphone.

- Sterilised by ethylene oxide (ETO).
- Clamp for greater safety in the sample sampling.
- Integrated part: the fins are injected on the needle, generating a single piece of greater safety.



Point in "Y" as secondary route for administration of medicinal products



*Bevelled needle, does not damage the silicone septum.
Made of stainless steel.*

Product/ code	Size (Ga) x length (mm)	Size (Ga) x length (inch)	Type/ configuration	Presentation
*1501-22-20-1	22 Ga x 20 mm	22 Ga x 3/4"	90° wings	10 Units
1501-22-25-1	22 Ga x 25 mm	22 Ga x 1"	90° wings	10 Units

** Check for availability.*