

Closed Urine Collection System



SILMAG

www.gruposilmag.com



BSI British Standards Institution
Quality Management System
N° FM 664794
ISO 13485

Closed Urine Collection System

Code 499



URINARY TRACT




This product has been designed for bladder emptying and sterile urine sample collection for diagnostic purposes in patients who require close urinary output follow-up.

- Silicone closure cap adaptable to the rigid connector.
- Clamp for burette emptying.
- Does not cause pain or friction injury.

Silicone
100%
Tube




100% silicone radiopaque tube with closed end.



The closed system reduces the risk of infection.

Graded burette for measuring urine output.



Clamp and rigid connector with plug for more safety in the socket shows.

Code	Size (FR)	Outer Ø (mm)	Inner Ø (mm)
499-04	4,0	1,40	0,75
499-06	6,0	2,00	1,00
499-08	8,0	2,70	1,50