

Percutaneous endoscopic gastrostomy kit (PEG)



SILMAG

www.gruposilmag.com



BSI British Standards Institution
Quality Management System
N° FM 664794
ISO 13485

PERCUTANEOUS ENDOSCOPIC GASTROSTOMY KIT (PEG)

Code 355

This product has been designed to administer food directly into the gastric cavity in patients unable to feed by mouth, patients with gastroesophageal reflux and upper gastrointestinal tract malformations, gastroesophageal surgery patients and temporary or permanent gastrostomy patients.



DIGESTIVE TRACK



Silicone
100%
Tube


Distal hemispherical button with free interior for feeding passage.

Dialing every 1 centimeter.



Greater loop opening.

- 100% silicone tube.
- Radiopaque head.
- Features a multifilament flexible loop for easy placement.
- Easy handling
- Less traumatic placement.
- Graduated rule for greater control in the guide insertion.
- Nitinol guide.
- Greater loop opening.




Kit:

- (A) 1 gastric tube.
- (B) 1 introducer needle.
- (C) 1 single and double top with cap.
- (D) 1 fixing ring with clip.
- (E) 1 scalpel.
- (F) 1 ergonomic pincers for clamping.
- (G) 1 sterile surgical field.
- (H) 1 flexible guide for placement

Code	Size (FR)	Outer Ø (mm)	Length (cm)
355-14-00	14,0	4,70	70
355-16-00	16,0	5,30	90
355-18-00	18,0	6,00	90
355-20-00	20,0	6,70	90
355-22-00	22,0	7,30	90
355-24-00	24,0	8,00	90



Ergonomic pincer and graduated rule for greater control in the insertion of guide.



100% silicone tube.