

Si/Flow

High flow cannula



www.gruposilmag.com



info@silmag.com.ar | (+54) 0358 4970205
Sarmiento 355 | Las Higueras | Córdoba | X5805BTG



Empresa fabricante de productos médicos, habilitada y certificada en Buenas Prácticas de fabricación de la Anmat, legajo N°229.



BSI British Standards Institution
Quality Management System
N° FM 664794
ISO 13485

Junior

Si/Flow

High flow cannula

Code 240



**AIRWAY
TUBES**



Si/Flow is designed for the flow requirements and anatomical characteristics of newborns, infants and children receiving high nasal flow therapy (HFN).

Silicone
100%
Cannula

Soft

Silicone 100% Product.
Rounded, atraumatic bevel.



Practical

Connections compatible with different brands on the market.
Rapid identification of measurements by color and iconography.

Comfortable

It provides freedom of movement. Its fixation system allows a total adaptation to the anatomy of the patient.



Si/Flow Junior



Code	Size	Patient (Kg)	Maximum flow (l/min)	Ø Outer nostrils (mm)	Length of cannula (cm)	Colour
240-XS	XS	0,5 – 2,5	5	2,3	25	Light blue
240-S	S	1 – 4	8	2,6	25	Red
240-M	M	1 - 10	10	3,0	30	Yellow
240-L	L	3 - 20	20	3,5	35	Violet
240-XL	XL	5 - 30	25	4,0	40	Green

